


## Electric Youth Bed

Model # 115E-H-A



### Features:

- Electric bed utilizes hand control for positioning
- Full safety length bedside rails that meet 2<sup>3</sup>/<sub>8</sub>" spacing requirements for preventing entrapment
- 5" braking casters
- Solid pan sleeping surface (non-helical design)
- Includes 4 IV hook holes and 4 urine bag hooks
- 5" Perfed Comfort Mattress

### Specifications

- Overall Dimensions: 72.5"L x 36"W
- Height Adjustment: 23" - 33"
- Sleeping Surface: 33.5" x 64.5"
- Weight Capacity: 350 lbs.

Phone | 716.759.7200  
Toll Free | 800.274.2742  
Fax | 716.759.0700

Email | [info@nkmedicalproducts.com](mailto:info@nkmedicalproducts.com)

[www.novummed.com](http://www.novummed.com) • [www.nkmedicalproducts.com](http://www.nkmedicalproducts.com)