



NOVUM

Medical Products

Manual Youth Bed

Model # 115M-H-A



Features:

- 3 Manual Cranks: Head Gatch, Knee Gatch, Hi-Lo Positioning
- Full safety length bedside rails that meet 2 $\frac{3}{8}$ " spacing requirements for preventing entrapment
- 5" braking casters
- Solid pan sleeping surface (non-helical design)
- Improved crank mechanism
- Includes 4 IV hook holes and 4 urine bag hooks
- 5" Perforated Comfort Mattress (NOT INCLUDED)

Specifications

- Overall Dimensions: 72.5"L x 36"W
- Height Adjustment: 23" - 33"
- Sleeping Surface: 33.5" x 64.5"
- Weight Capacity: 350 lbs.

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700

Email | info@nkmedicalproducts.com

www.novummed.com • www.nkmedicalproducts.com