



Millennium Bassinet

Model 9000X17T



Features:

- No Basket Required: Built in Break Resistant Panels
- Larger Sleeping Surface can accommodate twins
- All Stainless Steel Construction
- Easy to Operate Head Elevation Crank
- Removable Sleeping Surface for Easy Cleaning
- Deep 6" Drawer for Storage
- Identification Card Holder on End
- Lifetime Warranty on all steel components. 1 Year on all other parts.

Specifications:

- Overall: 43"H x 31"W x 17.5"D