

Maternity Bassinet

Model: 9201-C




Aesthetically pleasing while also meeting your practical needs



Features:

- Constructed of medical grade laminate
- Quiet simple reflux positioning
- Drawer Dimensions: 9”H x 28”W
- Includes two 10”x10” pull-out decks for extra work space as well as a side chart holder
- 5” Locking Casters
- 1 Year Warranty
- **Overall Dimensions**
38.5”H x 31.5”W x 18.5”D

Finish Options



Prestige Maple

Candlelight

Classic Cherry

Biltmore Cherry

Truffle