



Fetal Monitor Cart

Model 9202-C

Medical grade case carts constructed of aesthetically pleasing melamine material to create a warm medical environment. Designed to meet all of your practical needs with their various compartments and drawers




Features:

- Easy to Clean
- Extremely Durable
- Easy to Maneuver
- Custom Dimensions available
- 5" Casters, two locking and two non-locking
- Soft close drawer slides
- 2 Pull Out Tray, 3 Drawers
- Cord Cut-Out to allow for easy use as a monitor stand
- 5 Year Limited Warranty

Dimensions:

- Height: 43"
- Width of Base: 25"
- Width of Top: 23.8"
- Depth of Base: 25"
- Depth of Top: 21.8"
- Pull-Out Trays: 20"W x 14"D
- Drawer Faces: 21.8"W x 5.8"H

Finish Options:



Prestige Maple Candlelight Classic Cherry Biltmore Cherry Truffle