



NOVUM

Medical Products

Blanket Warmer

Model DC350



Designed for use in the operating room or labor and delivery rooms to keep patients warm throughout the entire procedure, reducing surgical site infection risks.

Features:

- Single-chamber blanket warming cabinet
- White epoxy coated steel exterior casing and interior insert
- Single pane, energy efficient, glass window
- Warming technology uses multiple zone heating to monitor surfaces ensuring safe temperature ranges
- Temperature Range: 90-160°F (Operates in Celsius or Fahrenheit)
- Low energy consumption and minimal heat loss due to balanced heating system
- Interior LED lighting with intensity and off settings, casts comforting blue glow when on
- Energy Efficient warming shut-off system that is separate from electronic control, prevents overheating
- Electronic control included
- Cabinet stores current programming so in case of power failure it will return to most current programming
- Aesthetically pleasing design with soft corners and curved accents



Specifications

Dimensions: 28"H x 18.5"W x 22"D

Capacity: 10-12 Blankets

Weight: 91 lbs

Inner capacity: 3.5 ft³

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com