



# NOVUM

## Medical Products

# Blanket Warmer

## Model DC400



*Designed for use in the operating room or labor and delivery rooms to keep patients warm throughout the entire procedure, reducing surgical site infection risks.*

## Features:

- Energy efficient WarmSafe Zone Heating monitors interior surface temperatures and balances heat throughout the cabinet with low energy usage and minimal heat loss.
- Silent operation with no fans or moving parts. No lint traps required.
- Convenient hands-free, push button door is designed to make loading easier.
- Single pane, energy-efficient, e-coated glass window in door.
- Door is fully gasketed and hinged on the right side of the unit.
- Non-skid rubber feet are standard, optional casters.
- Easy to clean epoxy-coated exterior and interior surfaces.
- Programmable control
- Interior lighting with a comforting blue glow with 2 intensity settings and off mode.
- Stackable design allows for customized configurations in a small footprint.
- Can be used on countertop or freestanding on floor.
- Adjustable temperature range of 32° - 71°C (90° - 160°F).
- Four digit L.E.D. display.
- On/off button & up and down adjustment buttons.
- Actual temperature button.
- Built-in speaker for audible feedback.
- Integrated button lock-out feature.



## Specifications

- Dimensions: 22"H x 24"W x 25"D
- Capacity: 14 - 18 blankets
- Weight: 85.5 lbs
- Inner Capacity: 4.0 ft<sup>3</sup>

**Safety:** In the event of a power failure the cabinet will remember its programming and begins to operate as before when power is restored.

**Safety:** A warming shut-off system, separate from the electronic control, prevents overheating.

Phone | 716.759.7200  
Toll Free | 800.274.2742  
Fax | 716.759.0700  
Email | [info@novummed.com](mailto:info@novummed.com)  
[www.novummed.com](http://www.novummed.com)