



## Blanket Warmer

Model DC750

### Features:

- Energy efficient WarmSafe Zone Heating monitors interior surface temperatures and balances heat throughout the cabinet with low energy usage and minimal heat loss.
- Programmable control. Silent operation with no fans or moving parts. No lint traps required.
- Single pane, energy-efficient, e-coated glass window in door.
- Door is fully gasketed and hinged on the right side of the unit.
- Non-skid rubber feet are standard, optional casters.
- Easy to clean epoxy-coated exterior and interior surfaces.
- Includes a heated center shelf to further optimize heat distribution throughout the cavity.
- Interior lighting with a comforting blue glow with 2 intensity settings and off mode.
- Stackable with DC400 blanket warmer
- Can be used on countertop or freestanding on floor.
- Adjustable temperature range of 32° - 71° C (90° - 160° F).
- Four digit L.E.D. display.
- On/off button & up and down adjustment buttons.
- Actual temperature button.
- Built-in speaker for audible feedback.
- Integrated button lock-out feature.
- Safety: In the event of a power failure the cabinet will remember its programming and begins to operate as before when power is restored.
- Safety: A warming shut-off system, separate from the electronic control, prevents overheating.



*Designed for use in the operating room or labor and delivery rooms to keep patients warm throughout the entire procedure, reducing surgical site infection risks.*



### Specifications

- Dimensions: 35.5”H x 24”W x 24.3”D
- Capacity: 20 - 25 blankets
- Weight: 134 lbs
- Inner Capacity: 7.5 ft<sup>3</sup>