



Diaper Scale

Model # NK2001

Constructed of easy to clean stainless steel, this unit provides a compact, accurate battery powered scale. Complete with an auto shut-off feature to maximize battery life.



Features:

- Compact battery operated
- Stainless steel construction
- Powered by internal rechargeable battery
- Auto shut off to increase battery life
- Large LCD read-out screen
- Weighs as fine as 0.1 gram

Specifications:

- Units: lbs, kg, g, oz, lb:oz
- Display: LCD read-out
- Dimensions: 8¾”L x 6½”W x 2”H
- Plate Dimensions: 6½” x 6½”
- Battery Life: approx. 10 hours per charge life
- Controls: Zero, On/Off, Units

