



## Blanket Warmer

### Model EC1540



*Designed for large capacity patient warming needs to keep patients warm throughout the entire procedure, reducing surgical site infection risks.*

## Features:

- The durable, reliable stainless steel blanket warmers quickly and efficiently heat blankets to precise temperatures using the digital controls.
- Mobile with heavy-duty casters.
- Fan-less technology is designed to provide quiet warmth for patients and the absence of lint-clogged fans make units safe and easy to maintain.
- Double-chamber warming cabinet with stainless steel exterior.
- The appliance is warmed using a low-heat-density electrothermal cable array. The electrothermal cable is positioned in the floor and two sides of the warming appliance, providing even heating of the interior chamber.
- Epoxy-coated blanket support assembly and center shelf.
- Programmable heat settings to save energy costs.
- Adjustable temperature range of 37° - 93°C (98° - 180°F).
- Built in timer allows for efficiency and maintains necessary temperature levels.
- Four digit L.E.D. display.
- On/off button & up and down adjustment buttons.
- Actual temperature button.
- Integrated button lock-out feature.
- Safety: In the event of a power failure the cabinet will remember its programming and begins to operate as before when power is restored.
- Safety: A warming shut-off system, separate from the electronic control, prevents overheating.



Shown with optional bumper

## Specifications

- Dimensions: 72.75"H x 24"W x 32.5"D
- Capacity: 40 - 50 blankets
- Weight: 374 lbs
- Inner Capacity: 15.4 ft<sup>3</sup>  
7.7 ft<sup>3</sup> per chamber