



Blanket Warmer

Model EC750



Designed for use in the operating room or labor and delivery rooms to keep patients warm throughout the entire procedure, reducing surgical site infection risks.

Features:

- Durable countertop blanket warmer
- Mobile with optional heavy-duty casters, allowing nurses and to maximize their time.
- Our fan-less technology is designed to provide quiet warmth for patients.
- The absence of lint-clogged fans make units safe and easy to maintain.
- Single-chamber warming cabinet with stainless steel exterior and white epoxy-coated interior insert.
- Door is fully gasketed with a magnetic closure.
- WarmSafe™ heating system incorporates a multiple zone warming technology that heats where and when it is needed. All chamber surface temperatures are monitored, providing an efficient balance of heat, low energy consumption and minimal heat loss.
- Includes a heated center shelf to further optimize heat distribution throughout the cavity.
- Four (4) non-skid rubber feet are standard.
- Stacks with EC750 (itself) blanket warmer
- Adjustable temperature range of 32° - 71°C (90° - 160°F).
- Four digit L.E.D. display.
- On/off button & up and down adjustment buttons.
- Actual temperature button.
- Built-in speaker for audible feedback.
- Integrated button lock-out feature.
- Safety: In the event of a power failure the cabinet will remember its programming and begins to operate as before when power is restored.
- Safety: A warming shut-off system, separate from the electronic control, prevents overheating.



Shown with optional casters

Specifications

- Dimensions: 35.5"H x 24"W x 25.75"D
- Capacity: 20 - 25 blankets
- Weight: 134 lbs
- Inner Capacity: 7.5 ft³