

# Bassinet with Bottom Shelf

*Model NB-SSxB*



*Pictured with optional  
polycarbonate basket and  
mattress pad*

## Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3" swivel, ball bearing, non-marking casters with threadguard.
- 2 braking, 2 steering casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA
- Options:
  - Polycarbonate basket: Model 69576x1-P
  - Bassinet Mattress: Model 7028x21
  - Extender shelf: Model 7028EX

## Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Warranty: Lifetime on all steel components.  
1 Year on all other parts.