

Bassinet with Bottom Shelf

Model NB-SSxB



Pictured with optional polycarbonate basket and mattress pad

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3" swivel, ball bearing, non-marking casters with threadguard.
- 2 braking, 2 steering casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA

Options:

- Polycarbonate basket: Model **69576x1-P**
- Bassinet Mattress: Model **7028x21**
- Optional Static-free casters are also available.

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Bottom Shelf stands approximately 10" inches off the ground.
- Warranty: Lifetime on all steel components.
- 1 Year on all other parts.