

Bassinet w. Drawer & Bottom Shelf

NB-SSxD



*Pictured with optional
polycarbonate basket and
mattress pad*

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
 - 3" swivel, ball bearing, non-marking casters with threadguard
 - 2 braking, 2 steering casters
 - Trendelenburg may be created from either end of the bassinet
 - Deep, secure basket well
 - Quiet-closing drawers designed to keep sleeping infants undisturbed.
 - Made in the USA
- Options:
- Polycarbonate basket: Model 69576x1-P
 - Bassinet Mattress: Model 7028x21
 - Extender shelf: Model 7028EX

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Drawer: 26.25"L x 15.25"W x 5.75"H
- Bottom Shelf sits approximately 9" off the ground (5" inches without casters)
- Warranty: Lifetime on all steel components, 1 year on all non-steel parts.