



Bassinet with Drawer

Model NB-SSxD



Pictured with optional polycarbonate basket and mattress pad

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 steering casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA
- **Options:**
 - Polycarbonate basket: *Model 69576x1-P*
 - Bassinet Mattress: *Model 7028x21*
 - Extender shelf: *Model 7028EX*

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Drawer: 26.25"L x 15.25"W x 5.75"H
- Warranty: Lifetime on all steel components. 1 Year on all other parts.