



## Bassinet with End Drawer

*Model NB-SSxED*



### Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- Drawer operates on nylon ball bearing roller slides
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 steering casters
- Trendelenburg may be created from either end of the bassinet
- Made in the USA

### Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Drawer: 22"L x 12.75"W x 5.75"H
- Shelf: 15.25"L x 12.60"W
- Warranty: Lifetime on all steel components, and 1 Year on all other parts.

### Options:

Polycarbonate basket: *Model 69576x1-P*  
Bassinet Mattress: *Model 7028x21*  
Extender shelf: *Model 7028EX*