

Bassinet w. Open Cabinet

NB-SSxOC



Pictured with optional polycarbonate basket and mattress pad

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 steering casters
- Trendelenburg may be created from either end of the bassinet
- Deep, secure basket well
- Open cabinets provide easy-to-access storage.
- Made in the USA
- Options:
 - Polycarbonate basket: Model 69576x1-P
 - Bassinet Mattress: Model 7028x21
 - Extender shelf: Model 7028EX

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Warranty: Lifetime on all steel components, 1 year on all non-steel parts.

