

Bassinet with Open Cabinet

Model NB-SSxOC



*Pictured with optional
acrylic plastic basket and
mattress pad*

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 steering casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA
- **Options:**
 - Polycarbonate basket: Model 69576x1-P
 - Bassinet Mattress: Model 7028x21
 - Extender shelf: Model 7028EX

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Cabinet: 27"L x 16"W x 16.25"H
- Warranty: Lifetime on all steel components.
1 Year on all other parts.