

Wood Front Bassinet

NB-WFxCC-3DL



*Pictured with optional
polycarbonate basket and
mattress pad*

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3DL (3 Dimension Laminate) Seamless Thermofoil wood front construction with 3 layers of material for durability
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 swivel casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA

Options:

- Polycarbonate basket: *Model 69576x1-P*
- Bassinet Mattress: *Model 7028x21*
- Extender shelf: *Model 7028EX*

Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
- Cabinet: 27"L x 16"W x 16.25"H
- Warranty: Lifetime on all steel components. 1 Year on all non-steel parts.