

Pediatrics

Specialty
Steel

Furnishings

NOVUM

Medical Products

Wood Front Bassinet

Model NB-WFxCC-3DL



*Pictured with optional
polycarbonate basket and
mattress pad*

Features:

- Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability
- 3DL (3 Dimension Laminate) Seamless Thermofoil wood front construction with 3 layers of material for durability
- 3" swivel, ball bearing, non-marking casters with threadguard
- 2 braking, 2 swivel casters
- Trendelenburg can be created from either end of the bassinet
- Deep, secure basket well
- Quiet self-closing cabinet
- Made in the USA
- **Options:**
 - Polycarbonate basket: *Model 69576x1-P*
 - Bassinet Mattress: *Model 7028x21*
 - Extender shelf: *Model 7028EX*


Specifications:

- Overall: 31"L x 17.5"W x 37.75"H
 - Cabinet: 27"L x 16"W x 16.25"H
- Warranty:** Lifetime on all steel components.
1 Year on all other parts.

Finish Options:



Fusion Maple




Mahogany




Crossfire Java




Studio Uptown



Majestic Walnut



Mocha Cherry



Amati Walnut

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com