

Diaper Scale

NK2002-2



Features:

- Stainless Steel casing, Stainless Steel platter, and ABS plastic base.
- IP68 Certified stainless-steel battery-powered scale.
- Scale provides 100 hours on a full charge. Wireless charging pad is included.
- Scale is completely sealed, so it's dust-tight and can withstand water and other fluids.
- Automatic shut-off timer for when the scale isn't being used.
- Can quickly switch between lb, oz, kg, and g.
- Rubber feet maintain scale's stability.
- One year warranty on scale parts.

Specifications:

- Overall Dimensions: 11"D x 9.5"W x 4.7"H
- Platter Dimensions: 7.25"D x 9.25"W
- Choose between 6 lb, 12 lb, and 30 lb weight capacity models.
- Display: 6 Digit, 1.2" backlit LCD display
- NTEP Certified
- IP68 Certified Environmental protection