

## iSeries Overbed Tables

NV-iOBT



### Features:

- Easy to clean with thermofoil or highly durable, 3D Laminate seamless top.
- Anodized aluminum column and heavy duty steel base with color options
- Low profile U-base for easy mobility to fit with beds, stretchers and recliners (1.75")
- Pneumatic lift mechanism to adjust height
- Finish options available for table tops
- Optional Ultra Low Base, 1.34" to fit under beds with 5th wheels and other furnishings
- Low profile base casters

### Ultra Low VS Standard Base

<u>iSERIES OVERBED TABLES</u>	<u>NV-iOBT</u>
STANDARD GREY BASE	NV-iOBT-S-SG
ULTRA-LOW GREY BASE	NV-iOBT-UL-SG
STANDARD TAN BASE	NV-iOBT-S-ST
ULTRA-LOW TAN BASE	NV-iOBT-UL-ST

TABLE TOP SHAPE AND COLOR WILL BE PROVIDED IN (PARENTHESES).

### Specifications

- Overall Dimensions: 17"W x 32"L
- Height Range: 28" - 44"
- Flip Side Table: 13.5"
- Weight Capacity: 100 lbs (evenly distributed)

<u>I SERIES OVERBED TABLES</u>	<u>NV-IOBT</u>
STANDARD GREY BASE	NV-IOBT-S-SG
ULTRA-LOW GREY BASE	NV-IOBT-UL-SG
STANDARD TAN BASE	NV-IOBT-S-ST
ULTRA-LOW TAN BASE	NV-IOBT-UL-ST

TABLE TOP SHAPE AND COLOR WILL BE PROVIDED IN (PARENTHESES).

## Table Top Design Options



Rectangle



Kidney Shape



Peanut Shape

## Ultra-Low Base


The lowest overbed table base in the industry!



## Add-on Options



Urinal Hook



Storage Baskets



Dual Vanity  
with Mirror