



## Stretchers

Model NV9000



### Features:

- Centrally controlled brake system
- Collapsible side rails
- Dual side foot pedals
- Retractable 5th wheel steering system
- X-Ray compatible backrest
- Shock Absorbing, non-marking bumpers
- Base storage area for oxygen E-cylinder tank
- Foley bag holder
- Removable IV Pole
- 3" Mattress included

### Specifications:

- Dimensions: 83"L X 33"W
- Weight Capacity: 550 lbs
- Patient Surface: 75.7"x 24.5"
- Height Range: 23"- 33.5"
- Caster Diameter: 8"
- Fifth Wheel Diameter: 4"
- **Positioning:**
  - Head: 84°
  - Knee: 30°
  - Trendelenburg: 16.5°
  - Reverse Trendelenburg: 16.8°