

RC Classic Medical Recliners

- Recliners available with Manual & Powered recline, as well as Standard & Bariatric sizes.
- Solid Steel frame, and heavy duty casters.
- Foot Tray, Side Tables, and Push Bar come standard (Does not apply to RC300s).
- Available with Heat & Massage features.



Classic Economy RC Recliners

RC300 | RC300-S

- 3" Casters rather than the usual 5".
- Available with or without Swing-Arms.
- 300 lb. weight capacity.



Classic RC Recliners

RC301 | RC301-S

- 5" fixed front Casters with accessible locks.
- Available with or without Swing-Arms.
- 300 lb. weight capacity.

Classic Bariatric RC Recliners

RC301BAR | RC301-S-BAR

- 5" fixed front Casters with accessible locks.
- Available with or without Swing-Arms.
- 500 lb. weight capacity.

	HEIGHT:	WIDTH:	SEAT DEPTH:	SEAT WIDTH:	SEAT HEIGHT:	ARM HEIGHT:	WEIGHT CAP.:	TRENDELENBURG DEPTH:
RC300 RC300-S	46.5"	32.75" (51.25" w/ TABLES UP)	19.5"	22.75"	20.5"	27.75"	300 LBS.	75"
RC301 RC301-S	48.5"	32.75" (51.25" w/ TABLES UP)	19.5"	22.75"	22.5"	29.75"	300 LBS.	75"
RC301BAR RC301-S-BAR	49.5"	36.5" (55.5" w/ TABLES UP)	20.5"	26.75"	25.5"	29.75"	500 LBS.	76"