

Medical Recliners



Model RC-T-ECC



Model RC-T-ECCB

Features


- Infinitely adjustable back
- Automatic extending leg rest
- 3 Position Lockbar
- Non-porous blow molded fold down side tray
- Trendelenburg release with handle accessible from either side.
- Built-in ergonomic headrest
- Elastic open seat base
- RC-T-ECCB Bariatric Capacity 450 lbs.

Optional Features:

- Pivot Tables
- Heat/Massage or Heat Only
- IV Pole
- Chair Headrest Cover
- Arm Rest Covers
- Side Cushion
- CAL 133

Specifications

- | RC-T-ECC | RC-T-ECCB |
|----------------------------------------------|-----------|
| ■ Weight Capacity: 350 lbs. or 450 lbs. | |
| ■ Height: 48.5 in. | |
| ■ Width: 33.5 in. or 37.5 | |
| ■ Maximum Width: 51.5 in. or 55.5 in. | |
| ■ Width with Tray Down: 33.5 in. or 37.5 in. | |
| ■ Seat Height: 21 in. | |
| ■ Seat Width: 21 in. or 25 in. | |
| ■ Seat Depth: 20 in. | |
| ■ Seat Back Height: 32.5 in. | |
| ■ Footprint Recline Length: 76.5 in. | |



Medical Recliners



Model RC-T-ECCS



Model RC-T-ECCSB

Features

- Infinitely adjustable back
- Swing Away Arms
- Automatic extending leg rest
- 3 Position Lockbar
- Non-porous blow molded fold down side tray
- Trendelenburg release with handle accessible from either side.
- Built-in ergonomic headrest
- Elastic open seat base
- RC-T-ECCSB Bariatric Capacity 450 lbs.

Optional Features:

- Pivot Tables
- Heat/Massage or Heat Only
- IV Pole
- Chair Headrest Cover
- Arm Rest Covers
- Side Cushion
- CAL 133

Specifications

- | RC-T-ECCS | RC-T-ECCSB |
|------------------------------------------|------------|
| ■ Weight Capacity: 350 lbs. or 450 lbs. | |
| ■ Height: 48.5 in. | |
| ■ Width: 34 in. or 38 | |
| ■ Maximum Width: 54 in. or 58 in. | |
| ■ Width with Tray Down: 34 in. or 38 in. | |
| ■ Seat Height: 21 in. | |
| ■ Seat Width: 21 in. or 25 in. | |
| ■ Seat Depth: 20 in. | |
| ■ Seat Back Height: 32.5 in. | |
| ■ Footprint Recline Length: 76.5 in. | |