



Medical Recliners

Model RC-T-PCRB

Features

- Sleek, easily cleaned Urethane Arms
- Fully Upholstered Sides
- 3 Position Lockbar
- Industrial elastic suspension system

Optional Features:

- Heat/Massage or Heat Only
- Wood Arms
- Side Tables

Specifications

- Height: 47 in.
- Maximum Height: 47 in.
- Width: 30 in.
- Upright Depth: 39 in.
- Seat Height: 20.5 in.
- Seat Width: 26 in.
- Seat Depth: 19.5 in.
- Seat Back Height: 29.5 in.
- Armrest Height: 28.25 in.
- Weight Capacity: 450 lbs.

