



# Medical Recliners

### RC 300 Medical Recliner

Economy unit arrives assembled with 3" European casters

### RC301 Medical Recliner/Transport

Intended for patient transport, this unit comes equipped with 5" caster, one side table, push bar and a foot rest

### RC301BAR Bariatric Medical Recliner/Transport

Designed for transporting the heavier patient, this unit is standard with 5" casters, one side table, push bar and foot rest.



Model RC301

*\*Electric Models Also Available with Same Specifications*



*Optional Infinite Recline Feature*

### Options:

- Infinite Recline
- IV Pole
- Side Table
- Swing Tray
- Breakaway Arm
- Drink Holder
- California Fire Code

	RC300	RC301	RC301BAR
Seat Height	22.5"	24"	24"
Seat Width	22.5"	22.5"	27.5"
Seat Depth	21"	21"	21"
Recline Sleep Length	74"	74"	74"
Overall Height	46"	48"	48"
Overall Width	32"	32"	37"
Depth Upright	41"	47"	47"
Chair Weight	220 lbs	220 lbs	220 lbs
Weight Capacity	300 lbs	350	500 lbs

### Fabric Color Options

