



Medical Recliner



Features:

- Available as a manual chair or electric powered
- Stationary arms
- 5" casters
- Wall-hugger recline mechanism
- Thermo-sensitive, pressure reducing foam
- Quick-release seat
- Durable reinforced steel, powder-coat frame
- Self-storing footrest

Specifications:

- Weight Capacity: 500 lbs.
- Back Height (from floor): 49"
- Overall Width (tables down): 37"
- Seat Width: 27.5"
- Steel Casters: 5"
- Trendelenburg Depth: 77"

Options

- Infinite Recline/Trendelenburg
- California Fire Code
- Side Table
- Power Chair Battery Backup
- Breakaway Arm
- Power Chair Dual Controllers
- Cupholder
- Heat or Heat and Massage

Fabric Color Options

