



Medical Cribs

Novum Hospital Cribs are designed to complement any medical setting, while ensuring the safety and comfort of the patient. We offer a selection of accessories and customizable options to meet your specific needs as well.

Standard Features:

- Sleeping surface options:
 - Gatch Deck, Trendelenburg Deck, Fowler Deck, Flat Pan Deck and Gatch/Trend Deck
- Safety extender included on all gatch decks
- Solid pan sleeping system in a non-helical design
- Easy access safety gates
- 5 Increments of Rail Stopping Points
- 5" Casters, 1 steering, 1 swivel and 2 braking
- Rail Spacing 2 3/8
- Crib Height: 62"
- Deck Height (from floor): 35"
- Available powder coated, chrome or stainless steel finish
- Lifetime Warranty
- Made in the USA



Optional Accessories:

- | | | |
|---------------------------|------------------------|-------------------------------------|
| ▪ Patient Weighing System | ▪ IV Pole with Bracket | ▪ Central Locking Casters |
| ▪ Cage or Plexi Top | ▪ Dialysis Gate | ▪ Undercrib Storage with O2 Cutouts |
| ▪ Crib Bumper Pads | ▪ Mattress | ▪ Head End Safety Gate |

Side-by-Side Specifications

INFANT SERIES			
Model #	Overall Size	Finish	Sleeping Surface
E1970CL	30"W x 44"L x 61"H	Epoxy	Flat Deck with Plexi End
C1970CL	30"W x 44"L x 61"H	Chrome	Flat Deck with Plexi End
E1970CLT	30"W x 44"L x 61"H	Epoxy	Trendelenburg Deck with Plexi End
C1970CLT	30"W x 44"L x 61"H	Chrome	Trendelenburg Deck with Plexi End
E1970FWCR	30"W x 44"L x 61"H	Epoxy	Crank Fowler with Plexi End
C1970FWCR	30"W x 44"L x 61"H	Chrome	Crank Fowler with Plexi End
E1980CL	30"W x 44"L x 61"H	Epoxy	Flat Deck
C1980CL	30"W x 44"L x 61"H	Chrome	Flat Deck
E1980CLT	30"W x 44"L x 61"H	Epoxy	Trendelenburg Deck
C1980CLT	30"W x 44"L x 61"H	Chrome	Trendelenburg Deck
E1980FWCR	30"W x 44"L x 61"H	Epoxy	Crank Fowler
C1980FWCR	30"W x 44"L x 61"H	Chrome	Crank Fowler
CHILD SERIES			
Model #	Overall Size	Finish	Sleeping Surface
E1981CL	30"W x 60"L x 61"H	Epoxy	Flat Deck
C1981CL	30"W x 60"L x 61"H	Chrome	Flat Deck
E1981CLT	30"W x 60"L x 61"H	Epoxy	Trendelenburg Deck
C1981CLT	30"W x 60"L x 61"H	Chrome	Trendelenburg Deck
E1981CG	30"W x 60"L x 61"H	Epoxy	Gatch Deck
C1981CG	30"W x 60"L x 61"H	Chrome	Gatch Deck
E1981CGT	30"W x 60"L x 61"H	Epoxy	Gatch/Trendelenburg Deck
C1981CGT	30"W x 60"L x 61"H	Chrome	Gatch/Trendelenburg Deck
E1981FWCR	30"W x 60"L x 61"H	Epoxy	Crank Fowler
C1981FWCR	30"W x 60"L x 61"H	Chrome	Crank Fowler
YOUTH SERIES			
Model #	Overall Size	Finish	Sleeping Surface
E1983CG	36"W x 72"L x 61"H	Epoxy	Gatch Deck
C1983CG	36"W x 72"L x 61"H	Chrome	Gatch Deck



NOVUM
Medical Products

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com