



## Visiting Nurse Scale

*Model # NK4001*



Constructed for durability, portability and accuracy. One-piece easy to use and sanitize design that is stable and reliable for any end use.

### Features:

- Designed for portability and accuracy
- One-piece sturdy and lightweight construction
- Portable carrying case available
- Reweigh Button: Zeroes the scale and reweighs the infant without having to move the baby

### Specifications:

- Max Weight: 45lb x 0.2oz
- Display: High LCD
- Accuracy: +/- 0.2oz
- Units: lb:oz, kg
- Battery Life: Four AA batteries (included)
- Overall Dimensions: 25"L x 13<sup>3</sup>/<sub>4</sub>"W x 3<sup>3</sup>/<sub>4</sub>"H
- Plate Dimensions: 25"L x 11"W