



## Visiting Nurse Scale

Model # NK4000



Adult Mode



Infant Mode

This portable scale is able to weigh infants, toddlers and adults all in one easy to use system. Its durable construction gives it the sturdiness to be ideal for visiting nurse applications as well as for use in doctors offices.

## Features:

- Portable, moisture resistant and light-weight unit
- Constructed of high-strength impact resistant plastic
- Moisture resistance allows for visiting nurse to carry through light rain or snow
- Design makes it easy to clean and maintain
- Two scales in one easy to manage system
- Battery powered
- Optional carrying case available

## Specifications:

- Max Weight: 400 lbs.
- Display: 4½ digit 7 segment LCD
- Accuracy: 0.2% +/- 1 digit of displayed resolution
- Units: pounds/ounces and kilograms
- Battery Life: est. 3,000 weigh cycles
- Adult Mode: 0.1 Lb (0.1 Kg)
- Infant Mode: ¼oz (0.01Kg)