



iSERIES

INNOVATE INTELLIGENT INSPIRED

iSeries Wallsaver Recliner

iRC2-S-LO20

Specifications:

- Tight seat and back / Removable upholstery
- 6" Wall proximity
- Upholstery Frames: Steel, Plywood & Hardwood
- Seat: High Density Polyurethane Foam Over Sheet Webbing
- Back: High Density Polyurethane Foam Over Sinuous Wire Springs
- Upholstered Covers: Attached w/ Velcro & Removable For Maintenance
- Antimicrobial Moisture Barrier Fully Encapsulates All Foam Areas
- Double Needle Stitching Standard

Dimensions:

- Width: 27"
- Depth: 34"
- Height: 41.5"
- Seat Height: 17.5"
- Arm Height 22.5"
- Seat Width: 21"
- Seat Depth: 21.5"

Other Specifications:

- Upholstery Frames: Steel, Plywood & Hardwood
- Seat: High Density Polyurethane Foam Over Sheet Webbing
- Back: High Density Polyurethane Foam Over Sinuous Wire Springs
- Upholstered Covers: Attached w/ Velcro / Removable For Maintenance
- Antimicrobial Moisture Barrier Fully Encapsulates All Foam Areas
- Double Needle Stitching Standard (Welting Optional - No Charge)



Chairs are pictured with optional arm caps

Chairs are available with casters