



# iSERIES

INNOVATE INTELLIGENT INSPIRED

## iSeries Wallsaver Recliner

*iRC2-S-WA20*



### Specifications:

- Tight seat and back / Removable upholstery
- 6" Wall proximity
- Upholstery Frames: Steel, Plywood & Hardwood
- Seat: High Density Polyurethane Foam Over Sheet Webbing
- Back: High Density Polyurethane Foam Over Sinuous Wire Springs
- Upholstered Covers: Attached w/ Velcro & Removable For Maintenance
- Antimicrobial Moisture Barrier Fully Encapsulates All Foam Areas
- Double Needle Stitching Standard
- Chair is Stationary unless optional casters are added.
- Optional Casters: 2 locking swivel, and 2 non locking, non swivel.

### Dimensions:

- Dimensions: 27"W x 34"D x 41.5"H
- Weight Capacity: 500 lbs.

**\*Chairs are pictured with optional arm caps\***

**\*Chairs are available with casters\***