

Glider Recliner

Model IRC-S-LO35



Features:

- Lever operated footrest
- Push on the arm operated back
- Chair automatically stops and locks gliding when reclined
- Seamless back
- Easy Clean out
- Contoured back with headrest and lumbar support
- Wide shoulders for comfortable patient positioning
- Stationary base
- Transport base available

Specifications:

- Dimensions: 46"H x 28.5"W x 35"D
- Seat: 19"H x 22"W x 21"D
- Fully Reclined Depth: 61"
- Weight Capacity: 350 lbs.
- Warranty: Limited 7 Year

Glider Recliner

Model IRC-S-LO35



Options:


Arm Caps



None



Wood
IRC-A-WA-A



Polyester
IRC-A-PA-A

Accessories



Push Bar
IRC-A-PB-A



Drainage
Bag Holder
IRC-A-DB-A



Restraint
Strap Loop
IRC-A-RS-A



IV Pole and
Holder
IRC-A-IV-A

*Rolling Foot
Tray
*Moisture
Barrier
* CAL-133

IRC-A-RF-A
IRC-A-MB-A
IRC-A-TB-A