



iSERIES Overbed Tables

Features:



- Easy to clean with thermofoil or highly durable, 3D Laminate seamless top.
- Anodized aluminum column and heavy duty steel base with color options
- Low profile U-base for easy mobility to fit with beds, stretchers and recliners (1.75)"
- Pneumatic lift mechanism to adjust height
- Finish options available for table tops
- Optional Ultra Low Base, 1.34" to fit under beds with 5th wheels and other furnishings
- Low profile base casters




Specifications

- Overall Dimensions: 17"W x 32"L
- Height Range: 28" - 44"
- Flip Side Table: 13.5"
- Weight Capacity: 100 lbs (evenly distributed)




Finish Options

Rigid Thermoformed/3DL




Thermoformed



Premium Upgraded Finish Options—100% Seamless RTF

Request a swatch kit for additional color matching quality.

Base Color Options



Gray



Cream

Table Top Design Options



Rigid Thermo/3DL



Thermoformed



Kidney Shape



Peanut Shape

Ultra-Low Base



The lowest overbed table base in the industry!



Add-on Options



Urinal Hook



Storage Baskets



Dual Vanity
with Mirror