

Overbed Tables

Model 120EV




Features:

- Raised urethane edge for spill containment
- Dual Cup Holders
- Flexible tabs grip the cup preventing spills
- Easily removed for cleaning
- Dual ergonomic handles
- Molded urethane provides comfortable grip
- Internal steel structure proves strength and durability
- Provides easy mobility
- Small storage area between the cup holders
- Base: Pneumatic lift mechanism provides height adjustment between 27" and 42"
- Internal bushing for smooth height adjustment
- Single handed height adjustment
- Available in both U and H style base design
- Column can be changed to mount on the right or left side
- Wide 28" opening allows easy access for most wheelchairs
- Opal powder coated base is standard
- Optional dual mirror vanity

Specifications:

- Top: 18" x 32"
- Height: 27" - 42"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year
- Ships fully assembled
- **120E same table as shown without vanity**



Finish Options



Bannister Oak

Clear Maple

English Oak

Figured Mahogany

Fusion Maple



Gunstock Walnut


Nepal Teak

Quartered Oak

Solar Oak

Wild Cherry

Standard Urethane Choices (Trim)



Sable

Putty

Opal

Black

Sandstone

Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

*Additional costs may apply to finish and trim upgrades