

## Single Top Overbed Tables

120E | 120EV




### Features:

- Raised urethane edge for spill containment.
- Available in alternate color choices.
- Optional vanity.
- Pneumatic lift mechanism provides infinite height adjustment between 27" and 42".
- Internal bushings for smooth, single-hand height adjustments.
- Column can be easily changed to mount on right or left side of base.
- Low height provides comfortable use from most patient room seating.
- Wide opening between base legs, allowing for easy access for most wheelchairs.
- Available in H and C style base designs.
- Opal powder coated base and column.
- Dual ergonomic handles are molded urethane to provide comfortable grip.
- Internal steel structure provides strength and durability.
- Provides easy mobility and helpful for caregivers, family, or patient.
- Cup holder's accommodates most commonly used beverage containers.
- Cup holder's flexible tabs grip the cup, preventing spills.
- Cup holders easily removed for cleaning.
- Storage area between cup holders perfect for pens, pencils, and other small objects.

### Specifications:

- Height: 27"-42"
- Top Surface: 18" x 32"
- Weight Capacity: 100 lbs.
- Weight: 55 lbs.
- Shipping Weight: Approx. 55 lbs.

***120EV same table as shown with vanity***



## Finish Options



Bannister  
Oak

Clear Maple

English Oak

Figured  
Mahogany

Fusion  
Maple



Gunstock  
Walnut

Nepal Teak

Quartered  
Oak

Solar Oak

Wild Cherry

### Standard Urethane Choices (Trim)



Black

### Standard Base Finish Choices




Brown



Gray



Opal



Black



Chrome

\*Additional costs may apply to finish and trim upgrades