

## Overbed Tables

### Model 120ESV




### Features:

- Dual Top/Split Top design
- Provides maximum space for food trays and personal items
- Upper includes a raised urethane edge for spill containment
- Lower top is hidden until easily activated with a single touch of a lever
- Tops slide and lock in either direction for added versatility and safety
- Dual cup holders
- Flexible tabs grip the cup preventing spills
- Easily removed for cleaning
- Dual ergonomic handles
- Molded urethane provides comfortable grip
- Internal steel structure provides strength and durability
- Provides easy mobility
- Small storage area between the cup holders
- Base: Pneumatic lift mechanism provides height adjustment between 29.5" and 44.5"
- Internal bushing for smooth height adjustment
- Available in both U and H style base design
- Column can be changed to mount on the right or left side
- Wide 28" opening allows easy access for most wheelchairs
- Opal powder coated base & column is standard

### Specifications:

- Upper top :18" x 32"
- Lower top: 17.5" x 26.5"
- Height: 29.5" and 44.5"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- 1 Year Warranty
- Ships fully assembled
- **120ES same table as shown without vanity**



## Finish Options



Bannister Oak

Clear Maple

English Oak

Figured Mahogany

Fusion Maple



Gunstock Walnut


Nepal Teak

Quartered Oak

Solar Oak

Wild Cherry

### Standard Urethane Choices (Trim)



Sable


Putty

Opal

Black

Sandstone

### Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome