

Overbed Tables

Model 120G



Features:

- High pressure laminate top
- Flush mounted vinyl t-mold
- Base: Gas cylinder lift mechanism provides quiet, smooth height adjustment between 27" and 42"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Available in both U and H style base design
- Column can be changed to mount on the right or left side
- Wide 28" opening allows easy access for most wheelchairs
- Opal powder coated base is standard
- Optional dual mirror vanity

Specifications:

- Top: 18" x 32"
- Height: 27" - 42"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 year
- ***120GV same table as shown with vanity***

Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut


Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Edge Choices



Brown

Opal

Black

Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

*Additional costs may apply to finish and trim upgrades