

Automatic Gas Lift Overbed Tables

120G | 120GV



Features:

Top:

- Walnut wood grain high pressure laminate on 3/4" particle board core.
- Flush mounted vinyl T-mold.

Height:

- Gas cylinder lift mechanism adjusts height between 30" - 45" smoothly and quietly.

Frame:

- Powder coated 11 and 12 gage welded tubular steel column and base.

Base:

- Heavy duty C-base with 28" clearance for use with chairs/wheelchairs.
- 1.5" x 2.5" tubular steel.
- Column can be mounted on left or right side.


Casters:

- 2" twin wheel

Specifications:

- Height: 30" - 45"
- Top Surface: 18" x 32"
- Weight Capacity: 100 lbs.
- Weight: 59 lbs.

120GV same table as shown with vanity



Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut

Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Urethane Choices (Trim)



Black

Standard Base Finish Choices



Brown



Gray



Opal



Black



Chrome

*Additional costs may apply to finish and trim upgrades