

Overbed Tables

Model 120GSV



Specifications:

- Upper top: 18" x 32"
- Lower top: 17" x 26"
- Height: 29.5" and 44.5"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year
- Ships fully assembled
- **120GS same table as shown without vanity**

Features:

- Dual Top/Split Top design
- Provides maximum space for food trays and personal items
- High pressure laminate top
- Resilient t-molding edge to resist scratching, chipping and stains
- Upper table slides either forward or back allowing the lower top to move into the opposite position and lock into place for added versatility and safety
- 2 slide activation tabs are located on column side of the table top
- Base: Gas cylinder lift mechanism provides quiet and smooth infinite height adjustment between 29.5" and 44.5"
Available in both U and H style base design
- Column can be changed to mount on the right or left side
- Wide 28" opening allows easy access for most wheelchairs
- Opal powder coated base is standard
- Optional dual mirror vanity

Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut


Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Edge Choices



Brown

Opal

Black

Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

*Additional costs may apply to finish and trim upgrades