

Automatic LTC Overbed Tables

121A & 123A




Features:

- Table top is walnut wood grain high pressure laminate on 3/4" particleboard core.
- Flush mounted T-mold edge.
- Chrome plated 12 and 14 gauge welded tubular steel column and base.
- 13 gauge stamped steel support bracket.
- Table height ranges from 29.5" - 44.5", adjustable in 15 one-inch increments.
- 1.5" twin wheel casters.

Specifications:


- Top: 15" x 30"
- Height: 29.5" - 44.5"
- 1.5" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year

*123A same table as shown
with vanity*



Finish Options

*Additional costs may apply to finish and trim upgrades



Bannister
Oak




Clear Maple




English Oak




Figured
Mahogany



Fusion
Maple




Gunstock
Walnut
(Standard)




Nepal Teak



Quartered
Oak



Solar Oak



Wild Cherry

Edge Choices



Brown
(Standard)



Opal




Black

Base Finish Choices



Brown



Gray



Opal



Black



Chrome
(Standard)