

## Overbed Tables

### Model 124




### Features:

- High pressure laminate top
- Flush mounted vinyl t-mold edge
- Base: Spring assisted height mechanism which locks into place in 1" increments between 30" and 45"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Available in both U and H style base design
- Column can be changed to mount on the right or left side
- Wide 28" opening allows easy access for most wheelchairs
- Standard Chrome plated steel column & base
- Standard Gunstock Walnut Finish
- Optional dual mirror vanity

### Specifications:

- Top: 18" x 32"
- Height: 27" - 42"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year
- ***124V same table as shown with vanity***



## Finish Options



Bannister  
Oak




Clear Maple



English Oak



Figured  
Mahogany



Fusion  
Maple




Gunstock  
Walnut



Nepal Teak



Quartered  
Oak




Solar Oak



Wild Cherry

### Standard Edge Choices



Brown



Opal



Black

\*Additional costs may apply to finish and trim upgrades