

Heavy-Duty Auto Overbed Tables

124 & 124V




Features:

- Table top is walnut wood grain high pressure laminate on 3/4" particleboard core.
- Flush mounted T-mold edge.
- Chrome plated 12 and 14 gauge welded tubular steel column and base.
- 13 gauge stamped steel support bracket.
- Table height ranges from 30" - 45", adjustable in 15 one-inch increments.
- 2" twin wheel casters.

Specifications:


- Top: 18" x 32"
- Height: 30" - 45"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed

124V same table as shown with vanity



Finish Options

*Additional costs may apply to finish and trim upgrades



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut
(Standard)

Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Base Finish Choices



Brown

Gray

Opal

Black

Chrome
(Standard)