

Automatic LTC Overbed Tables

125 & 129




Features:

- Table top is walnut wood grain high pressure laminate on 3/4" particleboard core.
- Flush mounted T-mold edge.
- Chrome plated 16 and 17 gauge welded tubular steel column and base.
- 13 gage welded tubular steel column and base.
- Available with H or C base.
- Table height ranges from 28" - 45", infinitely adjustable.
- 1.5" twin wheel casters.

Specifications:

- Top: 15" x 30"
- Height: 28" - 45"
- 1.5" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year

***129 same table as shown
with vanity***



Finish Options



Bannister
Oak




Clear Maple



English Oak



Figured
Mahogany



Fusion
Maple




Gunstock
Walnut
(Standard)



Nepal Teak



Quartered
Oak



Solar Oak




Wild Cherry

Standard Base Finish Choices



Brown



Gray



Opal



Black



Chrome
(Standard)

*Additional costs may apply to finish and trim upgrades