

Deluxe LTC Overbed Tables

126 & 126V




Features:

- Attractive and durable veneer top
- 1/4" spill protection lip
- Bullnose edgemold
- Heavy duty, height adjustable column with a slim and maneuverable base
- Infinite height adjustment between 28"-45"
- H style base design
- Standard Opal powder coated 16 & 17 gage welded tubular steel column and base
- 13 gage stamped steel support bracket
- Standard Gunstock Walnut finish
- Optional dual mirror vanity


Specifications:

- Top: 15" x 30"
- Height: 28" - 45"
- 50mm Twin Casters
- Weight capacity: 100 lbs evenly distributed

***126V same table as shown
with vanity***



Finish Options



Natural Pear

Honduras
Mahogany

Fusion
Maple




Gunstock
Walnut
(Standard)

Gray Nebula

Quartered
Oak

Wild Cherry

Standard Base Finish Choices



Brown

Gray

Opal
(Standard)

Black

Chrome

*Additional costs may apply to finish and trim upgrades