

## Overbed Tables

### Model 126




### Features:


- Attractive and durable veneer top
- 1/4" spill protection lip
- Bullnose edgemold
- Heavy duty, height adjustable column with a slim and maneuverable base
- Base: Spring assisted lift mechanism
- Infinite height adjustment between 28" - 45"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- H or U style base design
- Standard Opal powder coated steel column and base
- Standard Gunstock Walnut finish
- Optional dual mirror vanity

### Specifications:

- Top: 15" x 30"
- Height: 28" - 45"
- 2" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: Limited Lifetime on Steel Components
- **126V same table as shown with vanity**



## Finish Options



Natural Pear

Honduras Mahogany

Fusion Maple



Gunstock Walnut

Gray Nebula

Quartered Oak

Wild Cherry

### Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

\*Additional costs may apply to finish and trim upgrades