

Overbed Tables

128 & 128V




Features:

- Melamine laminate top
- Base: Automatic spring assisted lift mechanism
- Infinite height adjustment between 19.75" - 28"
- Designed for low beds with 3.25" clearance
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Available in both C and H style base design 24" opening
- Standard Chrome plated steel column & base
- Standard Gunstock Walnut
- Optional dual mirror vanity

Specifications:

- Top: 15" x 30"
- Height: 19.75" - 28"
- 40mm Dual Wheel Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year

***128V same table as shown
with vanity***



Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut

Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

*Additional costs may apply to finish and trim upgrades