

## Overbed Tables

### Model 132




### Features:

- Melamine laminate top
- 33 degree tilt in either direction
- Raised t-molded edge to prevent papers from falling when tilted
- Fingertip convenience for both tilt and height adjustment
- Base: Spring assisted lift mechanism
- Infinite height adjustment between 26.5"-39"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Standard Chrome finish steel column & base
- Standard Gunstock Walnut finish

### Specifications:

- Top: 15" x 30"
- Height: 26.5"-39"
- 1.5" Casters
- Weight capacity: 50 lbs evenly distributed
- Warranty: 1 Year



## Finish Options



Bannister  
Oak

Clear Maple

English Oak

Figured  
Mahogany

Fusion  
Maple



Gunstock  
Walnut

Nepal Teak

Quartered  
Oak

Solar Oak

Wild Cherry

\*Additional costs may apply to finish and trim upgrades