

Overbed Tables

Model 132




Features:


- Walnut wood grain melamine laminate top on 3/4" particleboard core.
- 33 degree tilt in either direction
- Raised t-molded edge to prevent papers from falling when tilted.
- Fingertip convenience for both tilt and height adjustment.
- Base: Spring assisted lift mechanism
- Infinite height adjustment between 26.5"-39"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- Standard Chrome finish steel column & base
- Standard Gunstock Walnut finish

Specifications:

- Top: 15" x 30"
- Height: 26.5"-39"
- 1.5" Casters
- Weight capacity: 50 lbs evenly distributed
- Warranty: 1 Year



Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut
(Standard)

Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

*Additional costs may apply to finish and trim upgrades